

**Percutaneous Sheath Introducer Kit with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters**

One: Sheath: 8.5 Fr. x 4" (10 cm) Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/Side Port and Tissue Dilator

One: Spring-Wire Guide: .035" (.89 mm) dia. x 17-13/16" (45 cm) (Straight Stiff Tip on One End - "J" Tip on Other)

One: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW

One: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 3 cc Luer-Lock Syringe

One: Catheter: 18 Ga. x 1-9/16" (3.97 cm) Radiopaque over 20 Ga. RW Introducer Needle with 3 cc Luer-Lock Syringe

One: 3 Swabstick Foil Pack 10% Povidone-Iodine

One: 5 ml Ampule HCl, 1% Lidocaine Solution

One: Cath-Gard® Catheter Contamination Shield: 80 cm

Three: Drape: 26" x 18" with adhesive

One: CSR Wrap

Five: Gauze Pad(s): 4" x 4"

One: Obturator: 8 Fr.

One: Scalpel: #11

One: Tape: Sterile

One: Stopcock: 3-way High-Flow

One: Suture: 3-0 Silk with Straight Needle One: SharpsAway® Disposal Cup